

CARDIOVIT AT-1 G2

The power of simplicity



SCHILLER
The Art of Diagnostics

CARDIOVIT AT-1 G2

- ∴ Easy patient data entry
- ∴ Outstanding signal quality for adult and paediatric ECG
- ∴ Lead reversal detection



WHEN TRADITION MEETS MODERNITY

Building on the solidly proven AT-1 and enhanced with the latest technology, the AT-1 G2 distinguishes itself by its outstanding signal quality and the newest interpretation algorithm.

User-friendliness is guaranteed with step-by-step workflow and easy patient data entry.

A colour screen and an easy-to-clean keyboard complete this intuitive and reliable electrocardiograph.

STATE OF THE ART

Thanks to a very high sampling frequency and a large frequency band, the CARDIOVIT AT-1 G2 provides optimal signal quality for adult and paediatric ECG.

The high-resolution thermal printer is compatible with roll and z-folded paper, accommodating customer preference. For data archiving PDF export on USB stick is also available.

ERGONOMIC DESIGN

The CARDIOVIT AT-1 G2 has been designed with ease-of-use in mind: patient data can be entered at the touch of a button; with a 5" colour display for ECG preview and simplified keyboard, taking an ECG has never been simpler.

STEP-BY-STEP WORKFLOW

1. Patient data entry (at the touch of a button)
2. Real-time signal preview and quality check, with lead reversal detection
3. Post-acquisition review

ADVANCED CLINICAL ALGORITHM

ETM: One of the industry's most renowned 12-lead ECG interpretation programs.

CCAA: SCHILLER's Culprit Coronary Artery Algorithm is used outside the hospital whenever a fast decision on chest pain is required immediately. It locates the obstruction site in the coronary artery and thereby the size of the cardiac area at risk.

1 The high-resolution thermal printer is compatible with roll and z-folded paper

2 Patient data entry with one touch of a button

3 ECG overview of 3 channels including colour-coded lead quality check

4 Lead reversal detection



EASY 1-2-3 STEPS

- 1: Enter patient data
- 2: Control signal quality and perform ECG
- 3: Review ECG before printing



Americas
SCHILLER Americas Inc.
Doral, Florida 33172

North America:
Phone +1 786 845 06 20
Fax +1 786 845 06 02
sales@schilleramerica.com
www.schilleramerica.com

Latin America & Caribbean
Phone +1 305 591 11 21
Fax +1 786 845 06 02
sales@schillerla.ch
www.schillerla.ch



Asia
SCHILLER Asia-Pacific / Malaysia
52200 Kuala Lumpur, Malaysia
Phone +603 6272 3033
Fax +603 6272 2030
sales@schiller.com.my
www.schiller-asia.com



Austria
SCHILLER Handelsgesellschaft m.b.H.
A-4040 Linz
Phone +43 732 70 99
Fax +43 732 757 000
sales@schiller.at
www.schiller.at



China
Alfred Schiller (Beijing) Medical Equipment
Co.,Ltd.
100015 Beijing, China
Phone +86 010 52007020
Fax +86 010 52007020-8016
info@schillerchina.com
www.schillermedical.cn



Croatia
SCHILLER d.o.o.
10000 Zagreb
Phone +385 1 309 66 59
Fax +385 1 309 66 60
info@schillerzg.hr
www.schiller.hr



France
SCHILLER Médical
F-67160 Wissembourg
Phone +33 3 88 63 36 00
Fax +33 3 88 63 36 49
info@schiller.fr
www.schiller-medical.com



France (distribution France)
SCHILLER France S.A.S.
F-77608 Bussy St Georges
Phone +33 1 64 66 50 00
Fax +33 1 64 66 50 10
contact@schillerfrance.fr
www.schiller-france.com



Germany
SCHILLER Medizintechnik GmbH
D-85622 Feldkirchen b. München
Phone +49 89 62 99 81 0
Fax +49 89 609 50 90
info@schillermed.de
www.schillermed.de



Hungary
SCHILLER Diamed Ltd.
H-1141 Budapest
Phone +36 1 383 4780 / +36 1 460 9491
Fax +36 1 383 4778
info@schillerhungary.hu
www.schillerhungary.hu



India
SCHILLER Healthcare India Pvt. Ltd.
Mumbai - 400 059, India
Phone +91 22 6152 3333 / 2920 9141
Fax +91 22 2920 9142
sales@schillerindia.com
www.schillerindia.com



Poland
SCHILLER Poland Sp. z o.o.
PL-02-729 Warszawa
Phone +48 22 843 20 89 / +48 22 647 35 90
Fax +48 22 843 20 89
schiller@schiller.pl
www.schiller.pl



Russia & C.I.S.
AO SCHILLER.RU
125040 Moscow, Russia
Phone +7 (495) 970 11 33
Fax +7 (495) 956 29 10
mail@schiller.ru
www.schiller.ru



Serbia
SCHILLER d.o.o.
11010 Beograd
Phone +381 11 39 79 508
Fax +381 11 39 79 518
info@schiller.rs
www.schiller.rs



Slovenia
SCHILLER d.o.o.
2310 Slovenska Bistrica
Phone +386 2 843 00 56
Fax +386 2 843 00 57
info@schiller.si
www.schiller.si



Spain
SCHILLER ESPAÑA, S.A.
ES-28232 Las Rozas/Madrid
Phone +34 91 713 01 76
Fax +34 91 355 79 33
schiller@schiller.es
www.schiller.es



Switzerland
SCHILLER-Reomed AG
CH-8953 Dietikon
Phone +41 44 744 30 00
Fax +41 44 740 37 10
sales@schiller-reomed.ch
www.schiller-reomed.ch



Turkey
SCHILLER TÜRKİYE
Okmeydanı-Sisli - Istanbul
Phone +90 212 210 8681 (pbx)
Fax +90 212 210 8684
info@schiller.com.tr
www.schiller-turkiye.com



United Kingdom
SCHILLER UK
Bellshill, ML4 3PR
Phone +44 1698 744 505
Fax +44 1698 744 474
sales@schilleruk.com
www.schilleruk.com

